




Racks & Cryo-boxes



Store your freezer storage boxes within organized freezer racks. METHER offers stainless steel racks to optimize your freezer space and makes it easy to find the items you desire. We offer variety of racks for our freezers. These racks are made of high-quality stainless steel and help you organize the samples inside your freezer. Our range consists of racks for both chest freezers and upright freezers.



Upright drawer rack

Strong stainless steel racks with drawers for easy and practical access to your boxes. The drawers slide smoothly and easily, which facilitates easy removal and insertion.

The upright drawer rack is equipped with a smart grip, which ensures easy handling and finding of necessary samples. When using racks with drawers, only the relevant drawer with the freezer boxes of interest is pulled-out.

All other boxes stay safely in the cold interior of the ULT freezer. The METHER racks are designed for smooth and ergonomic hand-lining, as well as to give a great overview of samples.



Model No.	N3J-254D	N3J-264D	N3J-255D	N3J-244D	N3J-243D	N3J-262D
Material	Stainless steel	Stainless steel	Stainless steel	Stainless steel	Stainless steel	Stainless steel
Box type	2 inch 9*9/10*10	2 inch 9*9/10*10	2 inch 9*9/10*10	2 inch 9*9/10*10	2 inch 9*9/10*10	2 inch 9*9/10*10
Box layout	5 x 4	6 x 4	5 x 5	4 x 4	4 x 3	6 x 2
Box Qty	20	24	25	16	12	12
Dimensions W*D*H(mm)	140*559*291	140*559*350	140*686*291	140*559*236	140*424*236	140*281*350

Chest front rack



The chest front rack is designed for optimal suiting of chest freezers. Due to the design and choice of material, these unique racks are extremely solid and long-lasting stainless steel. The box is removed from the front, thus ensuring a good overview. The handle in the top meanwhile ensures easy handling of the racks.

	N3J-2B1U	N3J-261U
Box layout	11 x 1	6 x 1



Cryo-boxes

2 inch box
133*133*53mm
9*9/10*10 wells

Paper cryo-box



PC cryo-box

PP cryo-box